Panasonic®

Operating Instructions (Household) The RF Ultrasonic Device Model No. EH-XR10

English.....	3
中文	23

Operating Instructions
(Household) The RF Ultrasonic Device
Model No. EH-XR10

Contents

Safety precautions.....	4	After use	17
Before use	9	Frequently asked questions	18
Parts identification	10	Separately sold Gel	19
Modes and their applications	11	Troubleshooting	20
Preparation.....	11	Battery life	22
How to use	13	Removing the built-in rechargeable battery.....	22
How to move the head in each mode	16	Specifications	22

Thank you for purchasing this Panasonic product.

Before operating this unit, please read these instructions completely and save them for future use.

Safety precautions

To reduce the risk of injury, loss of life, electric shock, fire, and damage to property, always observe the following safety precautions.

Explanation of symbols

The following symbols are used to classify and describe the level of hazard, injury, and property damage caused when the denotation is disregarded and improper use is performed.

DANGER

Denotes a potential hazard that will result in serious injury or death.

WARNING

Denotes a potential hazard that could result in serious injury or death.

CAUTION

Denotes a hazard that could result in minor injury or property damage.

The following symbols are used to classify and describe the type of instructions to be observed.

This symbol is used to alert users to a specific operating procedure that must not be performed.

This symbol is used to alert users to a specific operating procedure that must be followed in order to operate the unit safely.

DANGER

► Preventing accidents

Never use this appliance near electronic medical devices or home electrical therapy devices such as the following:

- Pacemakers and other electronic medical devices implanted in the body
 - Heart-lung machines and other life-supporting electronic medical devices
 - Wearable electronic medical devices, such as electrocardiographs and intravenous drip equipment
- Doing so may cause incorrect operation of electronic medical devices and home electrical therapy devices, resulting in accidents or causing patients to feel unwell.

WARNING

► To prevent skin or health problems

The following persons should not use this appliance:

- Individuals with heart disease
- Individuals with a nervous disorder in their face
- Individuals with hemorrhagic disease, hemophilia or other blood disorder
-  · Individuals undergoing treatment on their teeth
- Babies, toddlers or other child in the growth stage, or other individual who cannot indicate intent
- Individuals with aftereffects in the head due to a traffic accident, etc.
- Doing so may cause accidents or causing patients to feel unwell.

The following persons should not use this appliance:

-  · Persons with allergies
- Persons with sensitive skin
- Persons with skin diseases or atopic dermatitis
- Doing so may cause problems with the skin or body.

The following persons should not use this appliance:

-  · Pregnant women or within one month of childbirth.
- Doing so may cause skin trouble due to hormonal imbalance.

The following persons should not use this appliance:

-  · Individuals without sense of heat
- Doing so may cause burn.

Do not use this appliance in the following body area or following cases.

-  · When you do not feel well
- When you are taking medicine or when you are very tired
- Areas with skin damage, eczema or swelling
- When there are abnormalities on the skin such as dermatitis or extreme sunburns
- Body parts where epilation or shaving was performed on the same day
- Where surgery (plastic, cosmetic, etc.) has been performed
- Forehead (usable on forehead in RF mode only), lips, eyelids, ears, and neck
- Doing so may cause problems with the skin, body, or eyes.

WARNING

The following individuals must consult a doctor prior to use:

- Individuals with serious disease
- Individuals with acute disease
- Individuals with febrile disease
- Individuals with infectious disease
- Individuals with malignant tumors
- Individuals who, within the last 2 weeks, have undergone beauty treatment or some treatment by a clinic requiring a recovery period or have had treatment that caused severe redness in the skin
- Individuals with disturbance of perception due to a high degree of peripheral circulatory failure such as diabetes

If you feel pain or abnormality on your skin or other part of your body, immediately stop using this product.

- Failure to do so may cause accidents or skin/physical trouble.

The appliance has a heated surface. Persons insensitive to heat must be careful when using the appliance.

- Failure to do so may cause burn.

► **This product**

Do not use in proximity of flammable materials (such as alcohol, benzene, thinner, spray, nail polish remover).

- Doing so may cause fire or burn.

This product has a built-in rechargeable battery. Do not throw into fire or apply heat.

- **Do not charge, use, or leave in high temperatures.**
- Doing so may cause overheating, ignition, or explosion.

Do not store within the reach of children. Do not let them use it.

- Otherwise, small parts may be accidentally swallowed, thus causing accidents or injuries.

Do not modify or repair.

- Doing so may cause fire, electric shock, or injury.
- Contact Panasonic authorized service center for repair (battery change etc.).

Never disassemble except when disposing of the product.

- Doing so may cause fire, electric shock, or injury.

► **In case of an abnormality or malfunction**

Immediately stop using and remove the adaptor if there is an abnormality or malfunction.

- Failure to do so may cause fire, electric shock, or injury.

<Abnormality or malfunction cases>

- **The main unit, adaptor or cord is deformed or abnormally hot.**
- **The main unit, adaptor or cord smells of burning.**
- **There is an abnormal sound during the use of the main unit, adaptor or cord.**
- Immediately request inspection or repair at an authorized service center.

WARNING

► Power supply

Do not connect or disconnect the adaptor to a household outlet with a wet hand.

-Doing so may cause electric shock or injury.

Do not immerse the main body and AC adaptor in water or wash it with water.

Do not use or store the appliance in a bathroom or location with high humidity. Also, do not place in a location where water may spill (such as top of a vanity).

-Doing so may cause electric shock or fire.

Never use the appliance if the AC adaptor is damaged or if the power plug fits loosely in a household outlet.

Do not damage or modify, or forcefully bend, pull, or twist the cord or the power plug. Also, do not place anything heavy on or pinch the cord or the power plug.

-Doing so may cause electric shock or fire due to a short circuit.

Do not use anything other than the supplied AC adaptor for any purpose.

Also, do not use any other product with the supplied AC adaptor. (See page 10.)

-Doing so may cause burn or fire due to a short circuit.

Do not use in a way exceeding the rating of the household outlet or the wiring.

-Exceeding the rating by connecting too many plugs to one household outlet may cause fire due to overheating.

Always ensure the appliance is operated on an electric power source matched to the rated voltage indicated on the AC adaptor.

Fully insert the adaptor or appliance plug.

-Failure to do so may cause fire or electric shock.

Always unplug the adaptor from a household outlet when cleaning.

-Failure to do so may cause electric shock or injury.

Regularly clean the power plug and the appliance plug to prevent dust from accumulating.

-Failure to do so may cause fire due to insulation failure caused by humidity.

Disconnect the adaptor and wipe with a dry cloth.

► Take care of the following precautions as well

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.

Children should be supervised to ensure that they do not play with the appliance.

-Failure to do so may cause an accident or injury.

WARNING

-  **The supply cord cannot be replaced. If the cord is damaged, the AC adaptor should be scrapped.**
 - Failure to do so may cause an accident or injury.

CAUTION

► Power supply

-  **Do not allow metal objects or trash to adhere to the power plug or the appliance plug.**
 - Doing so may cause electric shock or fire due to a short circuit.

-  **Do not wrap the cord around the adaptor when storing.**
 - Doing so may cause the wire in the cord to break with the load, and may cause fire due to a short circuit.

-  **Disconnect the adaptor or the appliance plug by holding onto the adaptor or the appliance plug instead of the cord.**
 - Disconnecting by holding onto the cord may cause electric shock or injury.

-  **Disconnect the adaptor from the household outlet when not charging.**
 - Failure to do so may cause electric shock or fire due to electrical leakage resulting from insulation deterioration.

► Note the following precautions

-  **Do not allow the ultrasonic probe or RF poles to make contact with metals.**
 - Doing so may cause burn.

-  **Do not drop or hit against other objects.**
 - Doing so may cause injury.

-  **Do not use continuously for an extended period on the same area.**
 - Doing so may cause injury to the skin, infection, inflammation, bleeding, pigmentation or other trouble, or worsening of such symptoms.

-  **Before using, check that the ultrasonic probe is free from cracks, chips, deformities, and other damage.**
 - Failure to do so may cause injury to the skin.

-  **Excluding the arm used for a patch test, use only for care of the cheeks, face line, and forehead (in RF mode only).**
 - Failure to do so may cause skin damage, inflammation, etc.

Handling of the removed battery when disposing

DANGER

- The rechargeable battery is exclusively for use with this appliance. Do not use the battery with other products.
 - Do not charge the battery after it has been removed from the product.
 - Do not throw into fire or apply heat.
 - Do not hit, disassemble, modify, or puncture with a nail.
 - Do not let the positive and negative terminals of the battery get in contact with each other through metallic objects.
 - Do not carry or store the battery together with metallic jewellery such as necklaces and hairpins.
 - Do not use or leave the battery where it will be exposed to high temperatures, such as under direct sunlight or near other sources of heat.
 - Never peel off the tube.
- Doing so may cause overheating, ignition, or explosion.

WARNING

After removing the rechargeable battery, do not keep it within reach of children and infants.

- The battery will harm the body if accidentally swallowed.
- If this happens, consult a physician immediately.

If the battery fluid leaks out and comes in contact with the eyes, do not rub the eyes, and thoroughly rinse with clean water, such as tap water.

- Failure to do so may cause eye injury.
- Consult a physician immediately.

CAUTION

If the battery fluid leaks out and comes in contact with the skin or clothes, rinse off with clean water, such as tap water.

- Failure to do so may cause inflammation.

Before use

- The following symbol indicates that a specific detachable power supply unit is required for connecting the electrical appliance to the supply mains. The type reference of power supply unit is marked near the symbol.

► Usage frequency

- Use this appliance up to 12 minutes total for all modes (RF/Ultrasonic mode and RF mode) once during a day. Use up to 3 times per week and wait at least 1 day between each use.
- Overusing the appliance will actually hurt your skin.

► Skin precautions

- Before use, perform a skin allergy test using the conductive gel or some skin care product in place of the gel. (See page 11.)
- This quick test checks for any potential skin issues before you begin using the appliance normally.

- Please use this appliance with bare skin (after removing make-up and washing your face).
- Please check the condition of your skin while using this appliance.

When having abnormal sensation on skin, stop using the appliance immediately and consult a physician.

Skin issues can be attributable to stress, body condition, visceral disease, skin care products, and other factors. Seek professional medical attention if you have any skin issues.

▶ RF warmth level

- The ultrasonic probe increases in temperature as the RF warmth level is increased.

Select the desired RF warmth level according to your skin conditions and personal preferences.

When using the appliance for the first time or if any skin discomfort occurs, you are advised to adjust the RF warmth level to “LOW”.


▶ Sensors

- This appliance is equipped with a sensor to detect that the ultrasonic probe is making proper contact with the skin.
Over-care prevention sensor: Notifies the user when the ultrasonic probe has made contact with the skin for too long.
Fitted contact sensor: Notifies the user when the ultrasonic probe is not in contact with the skin.
 (See page 15.)

▶ Accessories

- Before using the appliance, be sure to remove any jewellery such as earrings on your face or around your neck.

Parts identification

A Main body

- ① RF Pole
- ② Head
- ③ Ultrasonic probe with warm sensation (Ultrasonic probe)
 - Provides ultrasonic energy and warmth.
- ④ Frame
- ⑤ RF & temperature indicator (Blue)
- ⑥ RF & temperature select switch
 - Works together with the RF level to adjust the warmth from the ultrasonic probe.
- ⑦ Mode indicator (Blue)

- ⑧ POWER/SELECT switch

- ⑨ Charging indication/ notification lamp (→) (Red)

- ⑩ Appliance socket

B AC adaptor (RE6-37) (The shape of the AC adaptor differs depending on the area.)

- ⑪ Adaptor
- ⑫ Power plug
- ⑬ Cord
- ⑭ Appliance plug

C Stand

Accessory

- ① Conductive gel (EH-4R01)
(See page 19.)

Modes and their applications

Select the mode appropriate to the desired application.

Mode	RF/Ultrasonic	RF
Area of application	 <p>Face line, cheeks, and mouth</p> <ul style="list-style-type: none"> • Do not use on the forehead. 	 <p>Forehead</p>
Function	RF × Ultrasonic energy × Vibration	RF
Usage time	10 min.	2 min.


Preparation

Skin (Allergy) test procedure

Make sure to check compatibility of this appliance with your skin as this product is normally used together with conductive gel. Skin care products you normally use can also be used in place of the gel.

This check must be performed to ensure you experience no skin issues as a result of use.

- Testing is required for each of your skin care products.

- 1 Clean the test area (inner arm) with soap and water, and dry it with a towel.**
- 2 Apply the conductive gel or skin care product.**
 - Apply a light amount just enough to cover the head.
- 3 Turn on the power and make contact with your skin.**
 1. Turn on the power and select the RF/Ultrasonic mode. (See page 13.)
 2. Once the appliance has completed the preheating process, press the ultrasonic probe onto your skin and slowly move the appliance in a circle for approximately 10 seconds. (See page 14.)
 - The head vibrates.

3. Remove the appliance from the skin and then press onto your skin and move the appliance for another 10 seconds.

- Once you have finished this procedure, allow the gel or skin care product to naturally dry.

4 Wait 48 hours without using the appliance.

- Do not allow the affected area to become wet and do not scrub this area.

If any skin issues occur within 48 hours

Do not continue to use the conductive gel or skin care product.


Charging

1 Insert the appliance plug into the appliance socket.

2 Mount the main body in the stand.

3 Insert the adaptor into a household outlet.

- Check that the charging indication/notification lamp (➡) is glowing.
- Charging is completed after approx. 2 hours.
The charging indication/notification lamp (➡) turns off.

4 After charging is completed, unplug the adaptor from the household outlet and unplug the appliance plug from the main body.

- (for safety and to reduce energy usage)
- Charging the battery for more than 2 hours will not affect battery performance.

<To check if charging is completed>

If you remove and reinsert the appliance plug while charging, the charging indication/notification lamp (➡) glows and turns off after approximately 5 seconds.

<Charging periods>

- Once the battery runs low, the charging indication/notification lamp (➡) will flash once per second. Once the battery is nearly depleted completely, an audible alarm will be triggered. Charge the appliance immediately.
- Once fully charged, the appliance can be used approximately 3 times. (in RF/Ultrasonic mode and RF mode)
However, the usage time may be reduced depending on the usage environment.

Notes

- When charging the appliance for the first time or when it has not been in use for more than 6 months, the charging time may change or the charging indication/notification lamp (➡) may not glow for a few minutes. It will eventually glow if kept connected.
- Recommended ambient temperature for charging is 15 °C – 35 °C. It may take longer to charge or the battery may fail to charge properly in extremely high or low temperatures. (See page 21.)

- If there is noise from the radio or other sources while using or charging the appliance, move to a different location to use the appliance.


How to use

Use the appliance on clean skin after thoroughly removing makeup and washing your face.

- A tone will be triggered after the POWER/SELECT switch is pressed or the sensor operates.
- Warming may fail when the ambient temperature is below 15 °C.

1 Press the switch and select a mode.


- Each press switches to the next mode. The selected mode flashes in blue.
- Select the mode appropriate to the desired application. (See page 11.)

- The appliance starts in RF/Ultrasonic mode each time the power is turned on.
- The appliance will not turn on when the AC adaptor is connected.
- Make sure to use the appliance with the frame attached.

2 Press the switch and select an RF warmth level to start the preheating process.

- Each press switches to the next level. (The selected mode flashes in blue.)
- The ultrasonic probe increases in temperature as the RF warmth level is increased.

Memory function

The previously selected RF warmth level is stored in memory so that the appliance can be used with the same setting.

This memory function will reset when the appliance has not been used for a long period of time or when the battery has ran out.

3 Apply the conductive gel while the appliance is preheating.

- We recommend using a conductive gel that helps transmit RF and ultrasonic energy.

1. Apply a macadamia nuts size (approximately 3 g) of gel into your hand.
2. Spread the gel around the areas of your face that you want to treat with this appliance as illustrated by the  areas in the figure.
 - Spread conductive gel gently on your face. Do not rub gel into skin.

Daily skin care products can also be used

- Use the conductive gel or skin care products to help transmit the RF and ultrasonic energy to the skin.
- Use skin care products that will provide a slippery surface for the appliance.

Skin care products with large amounts of oil can obstruct the transmission of RF and ultrasonic energy.


Gel, serum, and milky lotion Cream, toner, lotion, and oil

- If friction begins to develop during use, add some more skin care product to prevent any negative effects on the skin due to friction.

4 Preheating process completes (takes approximately 1 minute).

The preheating process will finish approximately 1 minute after turning on the power. A tone will be triggered, and indicators change from flashing to glowing states to indicate that the appliance is ready.


Notes

- The appliance can be used on the skin even while indicators are flashing. The lamps will begin glowing once contact is made with the skin.
- An alert tone will be triggered 15 seconds after the appliance is not in contact with the skin. Then, the power will automatically turn off after approximately 1 minute.

Changing the level after the preheating process completes

The level changes with each pressing of the .

- Select an appropriate level while paying attention to the condition of skin.


5 Once the preheating process completes, press the appliance onto your skin and begin moving it around. (See page 16.)

The power will be automatically turned off after each mode is completed.

- Press the POWER/SELECT switch to change the mode or turn off the power midway through.

How to apply the head properly

- Apply the appliance to the skin as both ultrasonic probe (a) and frame (b) are making complete contact.
- Press the appliance tightly to the skin as both ultrasonic probe and frame are contacting the skin in parallel.
The head vibrates except in RF mode.
- RF and Ultrasonic are only effective when the head is vibrating.

Note

- Head vibration produces sound, which is just normal operation.

Situations of an alert sound triggering followed by vibration stopping

Remove and press the ultrasonic probe on your skin again if either of the following notifications occur.

► Over-care prevention sensor

The ultrasonic probe has been in contact with your face for approximately 15 seconds.

► Fitted contact sensor

The ultrasonic probe has not made contact with the skin for approximately 15 seconds.

An alert will trigger and the power will turn off after approximately 1 minute if the appliance continues in either of these states.

How to move the head in each mode

RF/Ultrasonic mode - Approximately 10 minutes

(5 seconds per each movement shown by arrow in the figure)


English

Face line

Mouth to temples
(Approx. 8 min.)

Move the appliance as your cheeks are gently lifted upward.
(Repeat directions ① through ⑫.)


(A tone will be triggered 2 minutes before automatic power off.)

Mouth

(Approx. 2 min.)

Move the appliance upward around your face lines as your mouth is gently lifted upward.

(Repeat direction ⑬ through ⑱.)

(A tone will be triggered.)

Automatic power off


RF mode - Approximately 2 minutes

(5 seconds per each movement shown by arrow in the figure)

Forehead
(Approx. 2 min.)

Move from the center of forehead as indicated by ① and ② in the figure towards temples.

• Press the upper edge of the head above the eyebrows as indicated by ③, ④, ⑤, and ⑥ in the figure. Gently move the appliance upward.
(Repeat directions ① through ⑥.)

(A tone will be triggered.)

Automatic power off

Note

• If friction begins to develop during use, add some more conductive gel to prevent any negative effects on the skin due to friction.

After use


Cleanse the gel from the face.

- After cleansing off gel from the face, use your daily skin care products to condition your skin.
- If you used a daily skin care product instead of the conductive gel, you do not need to cleanse it off.

- After using the appliance, make sure to clean the main body and frame.
- Check if the power is turned off.

1 Pinch the frame to remove.

- If gel is making it difficult to remove the frame, wipe off excess gel and then remove the frame.

2 Clean the head.

Use a tissue or similar to wipe off gel or skin care product.


- Use a tissue, cotton swab, or similar to remove any gel or skin care product that has gotten into gaps.

3 Attach the frame.

Align the four notches (b) with the tabs (c) on the head and then press the frame onto the head until it clicks into place.

- Orient the frame so that the flat side (a) faces up.

English

4 Store the main body in the stand.

Notes

- The main body cannot be washed with water.
- If you are concerned that the main body is dirty, use a cloth soaked with soapy water to wipe the dirt. Wring out the cloth before wiping.
- Do not use alcohol, nail polish remover, or detergent (ex. liquid hand soap), etc. to clean the appliance. Doing so may cause failure, cracking or discoloration of parts.
- Make sure to remove all excess gel or skin care product. Failure to do so can cause hardening when this excess material dries.

Checking ultrasonic energy is transmitted in RF/ Ultrasonic mode

Transmission of ultrasonic energy cycles on and off for approximately 1 second for each phase, but this will not be felt by the skin. Use the following procedure to check transmission of ultrasonic energy.


1. Place a few drops of water on the ultrasonic probe (a).

- Orient the ultrasonic probe so that it is level.
- We recommend that you use water to check as it is more difficult if using gel.
- Do not allow any other parts beside the ultrasonic probe to come into contact with water.

2. Turn on the power and select the RF/Ultrasonic mode.

3. Use your finger to lightly press on the ultrasonic probe (a) until the ultrasonic probe is approximately the same height as the frame (b).

4. Confirm the transmission.

- The water droplets will vibrate if the appliance is operating properly. (The vibration cycle lasts for approximately 1 second.)

Frequently asked questions

Question	Answer
What is RF and ultrasonic energy?	<p>RF is an acronym for radio frequency. This is a type of high-frequency energy.</p> <p>▶ Ultrasonic energy includes sound frequencies over 20 000 Hz that cannot be heard by the human ear.</p>
Can I use this appliance on my hands and other areas of the body?	<p>▶ This product has been designed specifically for use on the face. Use this product only for beautifying the face.</p>
Can I use the appliance more frequently than once a day?	<p>Use this appliance up to 12 minutes total of all modes (RF/Ultrasonic mode and RF mode) once during a day.</p> <p>▶ Use up to 3 times per week and wait at least 1 day between each use. Overusing the appliance will actually hurt your skin.</p>
Why do I need to apply gel or skin care product?	<p>This material helps transmit the RF and ultrasonic energy to the skin.</p> <p>▶ (Any air between the ultrasonic probe and the skin causes the RF and ultrasonic not be transmitted.)</p>

Question	Answer
Can I use some other type of gel or milky lotion?	▶ While we recommend that use the conductive gel, serum, milky lotion, and gel can also be used. (See page 14.)
What type of skin care products should I use after using this appliance?	▶ Apply the skin care products you normally use on a daily basis after using this appliance.
Is it OK to charge the battery after every use?	▶ The battery in this appliance is a lithium-ion battery, and so charging it after every use will not affect battery service life.
What should I do if the gel causes staining or I develop skin issues?	Do not use the appliance until the skin condition clears. (We recommend that you seek professional medical attention.) ▶ You may have sensitive skin. Skin issues can also occur due to your physical condition, stress, body condition, and visceral disease.
What should I do if I feel tooth pain?	Use of this appliance may cause discomfort if you have any cavities. ▶ Do not use the appliance until you have had any cavities treated. (We recommend that you seek professional medical attention.)

Separately sold Gel

Conductive Gel (EH-4R01)

Troubleshooting

English

Problem	Possible cause	Action
The power does not turn on even after the POWER/SELECT switch is pressed.	<ul style="list-style-type: none"> ▶ The appliance has not been charged. ▶ The AC adaptor is connected. 	<ul style="list-style-type: none"> ▶ Charge the appliance. (See page 12.) ▶ Disconnect the AC adaptor before use.
The appliance can only be used for a few minutes after charging.	<ul style="list-style-type: none"> ▶ The battery has reached the end of its life. (approximately 3 years) (See page 22.) 	<ul style="list-style-type: none"> ▶ Please contact the store where purchased. (Customers cannot replace the battery as doing so may cause malfunction.)
The RF & temperature indicator, mode indicator, and charging indication/notification lamp (→) flashes once per second, and then the appliance stops operating.	<ul style="list-style-type: none"> ▶ The protection function in the main body has operated, and the appliance has stopped abnormally. 	<ul style="list-style-type: none"> ▶ Turn off the power and turn it on again after some time has passed. If the symptom occurs repeatedly, contact an authorized service center.
The RF & temperature indicator turns off during use.	<ul style="list-style-type: none"> ▶ The power will be automatically turned off after each mode is completed. 	
	<ul style="list-style-type: none"> ▶ After the preheating process completes, the power automatically turns off if no contact with the skin has been made for approximately 1 minute. 	
	<ul style="list-style-type: none"> ▶ A sensor has operated. 	<ul style="list-style-type: none"> ▶ Confirm that you are making proper contact with the skin. (See page 15.)
	<ul style="list-style-type: none"> ▶ There is no power remaining in the battery. 	<ul style="list-style-type: none"> ▶ Charge the appliance. (See page 12.)

Problem	Possible cause	Action
The charging indication/ notification lamp (➡) does not turn on or goes out while the appliance is charging.	<ul style="list-style-type: none"> ▶ This lamp turns off after charging is completed. (See page 12.) ▶ Make sure that the adaptor or appliance plug are inserted completely into place. 	
The charging indication/ notification lamp (➡) flashes twice per second while the appliance is charging.	<ul style="list-style-type: none"> ▶ This is an indication of abnormal charging. Make sure that you are charging the appliance in the recommended ambient temperature between 15 °C and 35 °C. 	<ul style="list-style-type: none"> ▶ If the symptom occurs repeatedly, contact an authorized service center.
The appliance is hot.	<ul style="list-style-type: none"> ▶ The appliance will be warm during use and charging. This is not indicative of any failure. 	
The appliance does not vibrate in RF mode.	<ul style="list-style-type: none"> ▶ The appliance does not vibrate in RF mode, but this is not a malfunction. 	

If the problems still cannot be solved, contact the store where you purchased the unit or a service center authorized by Panasonic for repair.

Battery life

Approximately 3 years.

If the operating time is significantly shorter even after a full charge, the battery has reached the end of its life. (The battery life may change due to usage and storage conditions.)

The battery in this appliance is not intended to be replaced by consumers. Have the battery replaced by an authorized service center.

Removing the built-in rechargeable battery

Remove the built-in rechargeable battery before disposing of the appliance.


The battery is to be disposed of safely.

Please make sure that the battery is disposed of at an officially designated location if there is one.

This figure must only be used when disposing of the appliance, and must not be used to repair it. If you dismantle the appliance yourself, it may cause malfunction.

- Remove the appliance from the AC adaptor.
- Press the POWER/SELECT switch to turn on the power and then keep the power on until the battery is completely discharged.

- Perform steps ① to ④ and lift the battery, and then remove it.
- Take care not to short circuit the positive and negative terminals of the removed battery, and insulate the terminals by applying tape to them.

For environmental protection and recycling of materials

This appliance contains a Li-ion battery.

Please make sure that the battery is disposed of at an officially assigned location, if there is one in your country.

Specifications

Power source	See the name plate on the AC adaptor. (Automatic voltage conversion)
Charging time	Approx. 2 hours

This product is intended for household use only.

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© Panasonic Corporation 2018-2020
Panasonic Corporation